

The widest Transitions® lens range available

Same great Transitions VI technology across all types of materials TL1.67, TL16, Airwear and Orma
Transitions lenses suit all wearers of single vision, progressives or bifocals

Transitions®
Healthy sight in every light™

PRODUCT	MATERIAL	TINT	SPHERE	ADDITION	CYL. TO
PROGRESSIVES					
Varilux Ipseo New Edition	TL1.67	Grey Et Brown	-12.00 to +9.00	+0.75 to +3.50	-6.00
Varilux Ipseo New Edition	Airwear	Grey	-10.00 to +6.00	+0.75 to +3.50	-6.00
NEW! Varilux Ipseo New Edition	Orma	Grey	-10.00 to +6.00	+0.75 to +3.50	-6.00
Varilux Physio Family	TL1.67	Grey Et Brown	-12.00 to +9.00	'+0.75 to +3.50 and +4.00 ¹	-6.00
Varilux Physio Family	TL16	Grey Et Brown	-12.00 to +8.00	'+0.75 to +3.50 and +4.00 ¹	-6.00
Varilux Physio Family	Airwear	Grey Et Brown	-10.00 to +6.00	'+0.75 to +3.50 and +4.00 ¹	-6.00
Varilux Physio Family	Orma 1.50	Grey Et Brown	-10.00 to +6.00	'+0.75 to +3.50 and +4.00 ¹	-6.00
Varilux Ellipse 360°	TL1.67	Grey Et Brown	-12.00 to +8.00	+0.75 to +3.00	-4.00
Varilux Ellipse 360°	TL16	Grey	-12.00 to +8.00	+0.75 to +3.00	-4.00
Varilux Ellipse 360°	Airwear	Grey Et Brown	-10.00 to +6.00	+0.75 to +3.00	-4.00
Varilux Ellipse 360°	Orma 1.50	Grey	-10.00 to +6.00	+0.75 to +3.00	-4.00
Varilux Ellipse	TL1.67	Grey Et Brown	-13.75 to +7.50	+0.75 to +3.00	-6.00
Varilux Ellipse	TL16	Grey	-12.00 to +8.00	+0.75 to +3.00	-6.00
Varilux Ellipse	Airwear	Grey Et Brown	-10.00 to +6.00	+0.75 to +3.00	-6.00
Varilux Ellipse	Orma 1.50	Grey	-10.00 to +6.00	+0.75 to +3.00	-6.00
Varilux Comfort DSx	TL1.67; TL16	Grey Et Brown	-12.00 to +8.00	+0.75 to +3.50	-6.00
Varilux Comfort DSx	Airwear	Grey Et Brown	-9.00 to +6.00	+0.75 to +3.00	-6.00
Varilux Comfort DSx	Orma 1.50	Grey Et Brown	-9.00 to +7.00	+0.75 to +3.50	-6.00
Varilux Comfort	TL1.67	Grey Et Brown	-13.00 to +7.50	+0.75 to +3.50	-6.00
Varilux Comfort	TL16	Grey Et Brown	-12.00 to +8.00	+0.75 to +3.50	-6.00
Varilux Comfort	Airwear	Grey Et Brown	-9.75 to +7.00	+0.75 to +3.00	-6.00
Varilux Comfort	Orma 1.50	Grey Et Brown	-9.00 to +7.50	+0.75 to +3.50	-6.00
Varilux Panamic	TL1.67	Grey	-11.50 to +7.50	+0.75 to +3.50	-6.00
Varilux Panamic	Airwear	Grey Et Brown	-10.00 to +7.00	+0.75 to +3.50	-6.00
Varilux Panamic	Orma 1.50	Grey Et Brown	-10.00 to +7.00	+0.75 to +3.50	-6.00
Definity and Definity Short	TL1.67	Grey Et Brown	-12.00 to +10.00	+1.00 to +3.50 and +4.00	-4.00
Definity and Definity Short	Poly	Grey Et Brown	-7.00 to +4.00	+1.00 to +3.00	-4.00
Definity and Definity Short	CR-39	Grey Et Brown	-10.00 to +6.00	+1.00 to +3.00	-4.00
Essilor Biofit FreeFrame and Essilor Biofit	TL1.67	Grey Et Brown	-12.00 to +8.00	+0.75 to +3.50	-6.00 ²
Essilor Biofit FreeFrame and Essilor Biofit	TL16	Grey Et Brown	-10.00 to +6.00	+0.75 to +3.50	-6.00 ²
Essilor Biofit FreeFrame and Essilor Biofit	Airwear	Grey Et Brown	-10.00 to +6.00	+0.75 to +3.50	-6.00 ²
Essilor Biofit FreeFrame and Essilor Biofit	Orma 1.50	Grey Et Brown	-9.00 to +6.00	+0.75 to +3.50	-6.00 ²
Essilor Ovation ADS	TL1.67	Grey Et Brown	-12.00 to +7.00	+0.75 to +3.50	-4.00
Essilor Ovation ADS	Airwear	Grey Et Brown	-10.00 to +6.00	+0.75 to +3.50	-4.00
Essilor Ovation ADS	Orma 1.50	Grey Et Brown	-8.00 to +6.00	+0.75 to +3.50	-4.00
Essilor Ovation	TL1.67	Grey Et Brown	-12.50 to +7.50	+0.75 to +3.50	-6.00
Essilor Ovation	Airwear	Grey Et Brown	- 6.25 to +7.50	+0.75 to +3.50	-6.00
Essilor Ovation	Orma 1.50	Grey Et Brown	-9.75 to +6.75	+0.75 to +3.50	-6.00
Essilor SmallFit	Airwear	Grey Et Brown	-10.00 to +6.00	+0.75 to +3.00	-6.00
Essilor SmallFit	Orma 1.50	Grey	-10.00 to +6.00	+0.75 to +3.00	-6.00
Essilor Natural	Airwear	Grey	-10.00 to +7.00	+0.75 to +3.00	-6.00
Essilor Natural	Orma 1.50	Grey	-8.00 to +7.25	+0.75 to +3.50	-6.00
Adaptar	Orma 1.50	Grey	-10.00 to +7.50	+0.75 to +3.50	-6.00
SEMI-FINISHED SINGLE VISION					
NEW! Essilor SV360°	TL1.67	Grey	-13.00 to +8.00		-4.00
NEW! Essilor SV360°	TL16	Grey	-13.00 to +8.00		-4.00
NEW! Essilor SV360°	Airwear	Grey Et Brown	-10.00 to +8.00		-4.00
AS TL 1.67	As TL1.67	Grey Et Brown	-17.00 to +8.00		-6.00
AS TL16	As TL16	Grey Et Brown	-13.00 to +10.00		-6.00
AS Airwear	As Airwear	Grey Et Brown	-12.00 to +9.00		-6.00
Airwear	Airwear	Grey Et Brown	-10.25 to +8.25		-6.00
Orma	Orma 1.50	Grey Et Brown	-10.50 to +7.00		-6.00
FINISHED SINGLE VISION LENSES					
NEW! FSV Crizal Forte	TL16	Grey	-2.00 to -6.00		-2.00
NEW! FSV Crizal Forte	As Airwear; Orma 1.50	Grey	-4.00 to +2.00		-2.00
FSV Orma Crizal Alizé	Orma 1.50	Grey	-4.00 to +2.00		-2.00
FSV Orma Crizal	Orma 1.50	Grey Et Brown	-4.00 to +2.00		-2.00
FSV AS Airwear HC	As Airwear	Grey	-6.00 to +4.00		-2.00
FSV Orma Supra	Orma 1.50	Grey Et Brown	-4.00 to +4.00		-2.00
BIFOCAL AND TRIFOCAL LENSES					
SF ST28	Orma 1.50	Grey Et Brown	-10.00 to +8.00	+1.00 to +3.00	-6.00
SF ST35	Orma 1.50	Grey Et Brown	-5.25 to +8.00	+1.00 to +3.00	-6.00
SF TRF 7x28 Orma	Orma 1.50	Grey Et Brown	-6.00 to +6.00	+1.50 to +3.00	-6.00

¹ Add. +0.75 to +3.00 for Varilux Physio Short 360°

² Cyl. -4.00 for FreeFrame




Essilor® Transitions® lenses with Crizal Forte™

Provide your patients with optimal clarity of vision outdoors and indoors while helping preserve the long term health of their eyes, so they will see better today and tomorrow.

PERFECT WITH
Crizal FORTE™





Transitions®




Why should you recommend Transitions® lenses with Crizal Forte™ to your patients?

- Transitions lenses with Crizal Forte adjust automatically to changing light for optimal clarity of vision outdoors and indoors
 - Reach dark state in approximately 1 minute (18% transmission and under, at 23°C)
 - Reach clear state in approximately 6 minutes (70% transmission and above, at 23°C)
- Transitions lenses with Crizal Forte enhance contrast and help reduce glare, relieving eye strain
- Transitions lenses with Crizal Forte block 100% of harmful UVA and UVB rays to help preserve the long term health of their eyes, so they will see better today and tomorrow

PERFECT WITH **CrizalFORTE™**

- Crizal Forte enhances the performance of Essilor Transitions lenses for better looking lenses
 - Best scratch resistance for long-lasting clarity of vision
 - Best antistatic and antismudge properties for cleaner lenses



True Light Control

Just as with a pair of clear lenses, your patient will need a second pair for specific uses such as driving or practicing sports in blinding glare condition.

For True Light Control, an additional pair of Essilor Polarizing lenses will provide both quality of vision and healthy sight in all types of light

High Resolution Vision



W.A.V.E. TECHNOLOGY™
WAVEFRONT ADVANCED VISION ENHANCEMENT

Transitions lenses are also available in High Resolution Vision for unrivalled contrast:

- Varilux Ipseo *New Edition*
- Varilux Physio f-360°
- Varilux Physio 360°
- Varilux Physio Short 360°
- Varilux Physio
- Varilux Ellipse 360°
- Essilor SV 360°



PROOF OF AUTHENTICITY!

Every pair of Essilor Transitions lenses is delivered with a Certificate of Authenticity (COA). Before you hand it to your patient, make sure you write down your Location identification Number (LIN) on the COA wrap. Every pair of lenses registered online by your patient earns you valuable points towards great rewards. Find out how by logging onto www.transitions.com/torp.

You can get your LIN number by contacting Transitions Optical Customer Service at 1-877-254-2590 or customerservicecanada@transitions.com.

TOTAL SATISFACTION
OR WE WILL GLADLY EXCHANGE THEM FOR OTHER ESSLOR LENSES*

CONDITIONS: the TOTAL SATISFACTION Program is applicable on any Essilor Transitions lenses coated with Crizal Forte D, Crizal Forte, Crizal EC, Trio EC, Trio or TD2 ordered from one of the Essilor Group laboratories and only in dissatisfaction cases. For the eye care professional, the exchange period is 3 months starting from the invoice date. The prescription (range, measurements) of the replacement lenses must be identical to the returned lenses. To obtain credit of the Essilor Transitions lenses, the eye care professional must provide to the laboratory: the Transitions lenses returned for dissatisfaction, a copy of the returned Transitions lenses invoice, and a copy of the replacement lenses invoice with the reason for return mentioned by the patient.

Visit transitions.com to learn more.